



SEPRAFLO[®] ION-EXCHANGE RESINS

Low-Cost Ion-Exchangers (Q & S) for Fast Protein Capture




- ✓ *50% Higher Flow than "Fast Flow" media*
- ✓ *Costs 60% less*
- ✓ *Excellent Binding capacity*
- ✓ *Non-compressible*
- ✓ *Cleanable in 1M NaOH*
- ✓ *1000 cycles*

Overview and Advantages

Septraflo[®] Resins are designed for process scale separation of biological molecules. They can be run at very high flow rates, typically up to 50% higher than commonly used resins from other suppliers without causing compression like agarose or cellulosic beads or fracture or breakage like silica or macroporous PSDVB based beads. They have large nominal pore size, up to 700 Å, that facilitates rapid mass transfer of large macromolecules and provides for high static and dynamic binding capacity. They are cleanable with 1N sodium hydroxide.

Septraflo[®] resins consist of macroporous beads with nominal 50 micron size. They are rigid and can withstand pressures in excess of 500 psi. They can be supplied by Septraflo in batches of 100 liters on up. The resins are made of a hydrophilic copolymer of glycidyl methacrylate which is used predominantly in making contact lenses. The Septraflo[®] resins are priced at a fraction of the price of competitive resins.

Ion-Exchange Binding Capacity of Septraflo[®]



The data above shows the binding capacity of Septraflo[®] resins when compared with other commercially available resins. The resins show comparable binding with a model protein BSA. Binding capacities will vary for different proteins based on their charge, size and molecular weight, PI, tertiary structure, and the pH and ionic strength of the solution and the competitive binding effect of other counterions in solution.



Application: Separation of Whey proteins using the Sepraflo® S Resin and Superflo® Column

The Process

The process for producing WPI, alpha-enriched WPI, beta-enriched WPI, lactoferrin, lactoperoxidase and glycomicropeptides is shown in the figure below:

In order to capture the proteins on the resin, the pH of the whey is adjusted to a low pH. This will give all the proteins a positive charge, which will bind it to the negatively charged Sepraflo®-S resin.

The pretreated and pH-adjusted whey is loaded on the column after the column has been equilibrated. Then the column is washed with water and the proteins are eluted by means of patented 'Opposing Salt and pH' gradients.

The process fractions are processed in an ultra- and diafiltration plant, where the protein is concentrated to 25-27% total solids, while the permeate containing the buffer salt solution is processed in a reverse osmosis plant, concentrating it to 2-3% solids for re-use.

The protein concentrate is dried in a spray dryer, which can be of the tall-form type using an inlet air temperature of 180°C.

The Sepraflo® resin combined with the Superflo® column process creates a cost effective next generation product opportunities, i.e. the isolation of Alpha Lactalbumin, Lactoferrin, Beta Lactoglobulin, GMP (GOS), IGs


Assumptions:

Volume: 250,000 liters of whey per day

Whey type: Sweet cheese whey, pre-treated, past., normal comp.

Processing time: 24 hours per day, 300 days a year


Ingredient/ Price	Traditional IEX	Sepraflo® Resins + Superflo® columns
Alpha lac \$20- 50/ Kg	Purity 20%- \$80% Cost \$18/ Kg	Purity 70% Cost \$10-14/ Kg
Beta Lac \$ 5/ Kg	Purity 40% Cost \$3/ Kg	Purity 90% Cost \$4/ Kg
Lactoferrin \$ 1000+/ Kg	Purity 90% Cost \$50/Kg	Purity 90%+ Cost \$30-\$50/ Kg






Application: Antibody Separation on Sepraflo® S Resin

The graph below shows an antibody separation on a Sepraflo® cation exchange resin. Elution is done with a salt gradient.



Flow Rate vs. Pressure Performance

The data below shows a comparison of the flow rates that can be achieved with various resins at varying pressures. Sepraflo® resins consistently show higher flow rates than all other resins.



Sepragen Has What It Takes

Sepragen's pioneering spirit of innovation brings to you the Sepraflo® Resins. You can now use these high performance resins to enable rapid, cost-effective purifications of antibodies and other proteins. Call 510-475-0650 to discuss your needs today!

Certificate of Analysis

DATE	December 9, 2016
PRODUCT	SEPRAFLO® S RESIN
LOT NUMBER	BGSS120901
QUANTITY	75 LITERS

PARAMETERS	SPECIFICATIONS	RESULTS	REFERENCE
Appearance	White Spherical Beads (Wet)	White Spherical Beads (Wet)	Visual Inspection
Total Exchange Capacity, Meq/ml (in H ⁺ form)	220-300 Meq/ml	282	QA/FPR/09
Moisture Content, % (in H ⁺ form)	62 - 70%	64.09	QA/FPR/05
Static Protein Capacity, mg/ml	65 - 85 mg/ml	78.90	QA/FPR-131
Particle Size Distribution:			
US Mesh			
+ 140 MESH	5% Maximum	4.70	QA/FPR/02
+ 270 MESH	85% Minimum	93.10	
- 270 MESH	10% Maximum	2.20	
Microscopic Analysis:			
Whole Beads	95% Minimum	99.00%	QA/FPR-06
Broken Beads	3% Maximum	1.00%	
Organic Leachability Analysis:			
DM Water Leachability, µg/ml	1 Maximum	Nil	FDA Publication 21 CFR 173.25
15% Ethanol Leachability, µg/ml	1 Maximum	0.40	
5% Acetic Acid Leachability, mg/ml	1 Maximum	0.82	

Specifications:

Matrix: Glycidyl Methacrylate Copolymer
 Functional Groups: Quaternary Amine and Sulphopropyl
 Pore Size: ~700 angstrom
 Particle Size: 20-100 microns (average 50 microns)
 Binding Capacity BSA : 70 mgs/ml



33470 WESTERN AVENUE, UNION CITY, CA 94587
 PHONE: 510-475-0650 FAX: 510-440-5395
 HTTP://WWW.SEPRAGEN.COM